



Nuova SB System S.r.l.

GANDY ETS

Endoscope Transit System



GANDY CARTS

Our transport carts can be used for both short-term storage and transport of endoscopes. They are characterized by different colors to signify clean, patient-ready scopes versus contaminated scopes that need reprocessing in order to PREVENT CROSS-CONTAMINATION.

Gandy ETS is designed for the storage & safe transport of seven endoscopes.



GANDY ETS

TECHNICAL DETAILS


CLEAN

Construction of cart	Stainless Steel 316
Cart Dimensions	145,4 cm X 74,9 cm X 53,1 cm (HxWxD)
Internal dimensions of Drawers	15,8 cm X 43,1 cm X 47,6 cm (HxWxD)
Color of drawers	Green
Weight	70kg




DIRTY


Construction of cart	Stainless Steel 316
Cart Dimensions	145,4 cm X 74,9 cm X 53,1 cm (HxWxD)
Internal dimensions of Drawers	15,8 cm X 43,1 cm X 47,6 cm (HxWxD)
Color of drawers	Red
Weight	70kg



GANDY ETS



Use with Basket



Use without Basket

- Gandy Endoscope Transit System (ETS) includes cart and drawers
- Usable also without additional baskets
- Made of stainless steel as an open design which is easy to clean and sterilize
- Compatible with stainless steel baskets developed and designed by Nuova SB System for “Gandy-Endoscope Treatment Products”
- Latex free
- Security lock for Endoscopes protection



Nuova SB System S.r.l.

Nuova S.B. System S.r.l.
Via Margherita de Vizzi Viganò 90 - 20092 - Cinisello Balsamo (MI)
Tel. 02.91474579 - Fax. 02.93664575
P.IVA e C.F 02223180395 - E-mail: info@sbsystem.it - PEC: nuovasb@registerpec.it