

Peerless™

HPLC Column Washing Pump



Features:

- Prompt and user friendly operation
- User friendly column Mounting
- Operation Modes: Continuous or timed
- All parameters and status simultaneously shown in single screen
- Simultaneous flushing of upto 5 columns

Your partner in chromatography...

- Automatic pump stoppage after completion of column flushing
- Inbuilt purge valve
- Inbuilt Manifold
- Automatic flow rate calculation for more than one column
- Inbuilt over current protection
- Negligible noise
- Compact table top assembly

Specifications:

- Flow range : 0.02ml/min – 9.99ml/min
- Pressure range : 0 psi – 6500 psi
- Automatic stop timer range : 1 min – 999 mins
- No. of column connections: 5(max)
- Max. current consumption : 3.5 A
- Minimum required bench space : 50cm X 25cm
- Required room temp range: 4 - 35° C
- Digital Pressure gauge : present
- Leak detection and high back pressure detection with alarm
- Auto flow stop mechanism in case of excess backpressure

Your partner in chromatography...

