



A99-1

Multi-functional
Electric Obstetric bed

SKF-C



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

DIMENSION PARAMETER

| | |
|-------------------|--------|
| External Length | 2000mm |
| External Width | 600mm |
| Table Height Min. | 650mm |
| Table Height Max. | 900mm |

TECHNICAL PARAMETERS

| | |
|-------------------------------------|--------------------|
| Reverse Trendelenburg/Trendelenburg | 10°/12° |
| Back Plate | 50° |
| Leg board Outwards | 90° |
| Main Voltage AC | 220V±22V, 50Hz±1Hz |
| Input power | 150±10% |
| Patient Weight Capacity | 185kg |

STANDARD ACCESSORIES

| | |
|-------------------------|------|
| Bed Frame | 1pc |
| Mattress | 1pc |
| Pull Handle | 2pcs |
| Armrest | 2pcs |
| Knee Curtch | 2pcs |
| Pedal Remote Controller | 1pc |
| Stainless Steel Basin | 1pc |
| Instrution Manual | 1pc |

QUALITY CONTROL



LEG BOARD OUTWARDS

The main bed and auxiliary bed are connected with each other, auxiliary table is easy for rotating and dismantling



PEDAL REMOTE CONTROLLER

Imported LINAK motor, Various actions are controlled by remote controller.




MATERIAL

The base seat and upright column covering are made of highquality 304 stainless steel.



STAINLESS STEEL BASIN

Equipped with stainless steel sewage basin and other various accessories



BRAKING SYSTEM

Foot treadle brake device which is easy for moving obstetric table.



POWDER SUPPLY

Internal battery maintaining power supply for 50-80 operations and with low battery warning function

PROFESSIONAL FUNCTIONS




RICH PARTS



Armrest



Pull Handle



Pedal Remote Controller



Knee Curutch



Stainless Steel Basin



INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.