



# A99-7


Multi-Functional  
Electric  
Obstetric Bed

SKF-C



**SAIKANG** / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE  
Stock Code: 870098



Height Adjustment

Back Plate up/down

Seat Plate up/down

Leg Plate Swing Out

## TECHNICAL PARAMETERS

External L*W*H	1980*720*(630-880)mm
Back Plate up/down	0-55° / 0-18°
Seat Plate up/down	0-35° / 0-5°
Leg Plate Swing Out	0-90°

## TECHNICAL CONFIGURATION

Bed Frame	1pc
Foot Pedal	1pc
Hand controller	1pc
Anesthesia Screen	1pc
Hand Bar	2pcs
Leg Support with Cushion	2pcs
Arm Support with Cushion	2pcs
Stainless Steel Sewage Basin	1pc
Caster	1pc

## ELECTRONICS

Voltage	100-240V
Frequency	50/60Hz
Plug	as required (hospital grade)

## QUALITY CONTROL



### MATERIAL

The base seat is plastic cover and height adjustment system is made of iron, environmental plastic can be sterilized in 100°C, tensile upto 30MPa

100c | 250MPa

### FUNCTIONS

The main bed and auxiliary bed are connected with each other, auxiliary table is easy for rotating and dismantling



### CONTROL SYSTEM

Handle controller which is easy for moving obstetric table.

### LEG SUPPORT

Double using leg support, can be used as boatstretcher



### MATTRESS

The mattress adopts high density memory sponge and is once formed, with no seam.



### FOOT STEP

Equipped with stainless steel sewage basin and other various accessories



### CONTROL SYSTEM

Good quality domestic motor, Various actions are controlled by remote controller.

## RICH PARTS

LEADING BRAND OF  
CHINA MEDICAL FURNITURE



Pull Handle



Armrest



Knee Curch



Pedal Remote Controller



Stainless Steel Basin





## INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

### JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901   
[www.chinasaikang.com](http://www.chinasaikang.com)

#### SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,  
Jiangsu Province,China.  
TEL: +86-512-58525850  
FAX: +86-512-58520713  
E-mail: [inquiry@jssaikang.com](mailto:inquiry@jssaikang.com)

#### SAIKANG (USA Branch)

Saikang Medical Technology, Inc  
17800 Castleton St Ste 180,  
City of Industry, Ca 91748  
Tel: +1 949648-2776

#### SAIKANG (UK Branch)

UK Healthward International Limited  
Owner(s) Address-DEPT 111,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +44 2032399738

#### SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.  
Viraj Complex Unit 3 Mombasa Road,  
PO Box 14262-00800, Nairobi.  
TEL: +254 733-704 690

#### SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC  
X22-S03 England Cluster,  
International City, Dubai UAE  
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.  
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.