



SKB038-2

PATIENT TROLLEY

Registration code: SKM-E





SAIKANG / Medical

Healthier
your life


SKB038-2

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External size (LxW)	1960x700 mm
▪ Table height adjustment	620-860 mm
▪ Back plate adjustment	0-75°



Technical configuration

<input checked="" type="checkbox"/> Crank system	1 set
<input checked="" type="checkbox"/> Back plate air spring	1 set
<input checked="" type="checkbox"/> Plastic side guardrail	1 set
<input checked="" type="checkbox"/> 5" covered castors	4 pcs
<input checked="" type="checkbox"/> IV pole	1 set
<input checked="" type="checkbox"/> IV pole prevision	4 pcs
<input checked="" type="checkbox"/> Mattress	1 pc

Guardrail

Virgin ABS material tuck away side rail with buffer. The guardrail is shaped on easy structure to hold and it can help the patient stand up easily by holding the rail.



Handle crank

Steel hanging crank, 100,000 times remaining intact, it can be fold under the bed when not in use.



Bed surface

The trolley surface is made of powder coated steel, smooth surface no welding scar.



Bed frame

The trolley base is made of powder coated steel.




Caster

5 inch TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts. Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500time.



Gas spring

Back rest lifting 0-75°, controlled by gas spring.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue