



NEED TO REPEATEDLY AND ACCURATELY DISPENSE POWDER?

- Consistent and Accurate
- Ergonomic and Easy to Use
- Ideal for Glove Boxes
- Typically $\leq 2\%$ RSD



VERSATILE

The dispensing volume for each handle is manually adjustable and multiple handle sizes can be used with each GamPette controller unit. This allows a wide powder volume range to be dispensed.



PRECISE AND REPRODUCIBLE

The relative standard deviation achieved on dispensed weights is dependent on powder properties, however, $\leq 2\%$ RSD is typically achieved.



ERGONOMIC

The lightweight handle allows easy and fatigue-free dispensing, ideal for glove box use. The handle is connected to the compact controller unit via a flexible vacuum pipe.



There's nothing else quite like Gamlen. Gamlen helps your formulation, production and manufacturing teams make better informed decisions every time.

MAKE IT BETTER. MAKE IT GAMLEN.