



## medifa 3000

Mobile functional furniture for OR,  
outpatient department and doctor's practice



we care.



Our mobile functional furniture provides safety and high quality in every respect. The range contains a variety of clever, mobile assistants for all functional areas of a clinic.

There is a focus on flexibility and hygiene in all of the products, whether they are furniture, seating or equipment. All of the furniture is made of stainless steel and thus corrosion-resistant, resistant to disinfectants, hygienic and easy to care for - throughout its entire service life.

**Your advantages at a glance**

- All products of the medifa assist series are made of ground stainless steel and have a very fine and hygienic surface.
- Manoeuvrable synthetic antistatic double castors ensure optimum mobility and are abrasion-resistant, non-discolouring and corrosion-resistant.
- The five armed footbase provides a high level of stability.
- Available in numerous equipment variants
- High level of stability.
- Welded without joints to eliminate corners and edges.
- The electrically conductive design of all relevant parts prevents static charging.
- Reliable service.

CONTENT	PAGE
Mayo instrument tables .....	4
OR swivel chairs and OR stools.....	6
Instrument trolleys .....	8
All-purpose and transport trolleys ...	10
Hook systems .....	12
Accessories for hook systems.....	15
Further OR furniture.....	21
Examination couches .....	24
Mobile stools.....	26
Upholstery colours .....	28
Further products for the OR .....	29



# Mayo instrument tables



## Hydraulic pump mayo instrument table, model Tübingen

**Art.-no. 30001248**

- Mobile 4-leg construction made of stainless steel, seamlessly welded
- Extra large table top made of stainless steel, with rounded corners and edges, 360° rotatable and lockable
- With built-in precision pump for height adjustment from 950 to 1350 mm
- 4 swivel castors Ø 125 mm, electrically conductive, 2 of them with brakes
- Max. load: 30 kg
- Dimensions: 1200 x 600 x 950 / 1350 mm
- Weight: 43 kg

## A handle on safety

It is very important to us to meet special requirements in the OR. The products of medifa assist therefore stand for maximum functionality and optimum hygiene conditions.



## Hydraulic pump mayo instrument table

**Art.-no. 30001246**

- Heavy, tilt resistant 3-foot construction made of stainless steel with built-in precision pump
- Stable, deep-drawn table top made of stainless steel, with rounded corners and edges, 360° rotatable and lockable, 640 x 440 mm
- Height adjustment by foot pedal from 950 to 1350 mm

- 3 double swivel castors Ø 75 mm, 1 of them electrically conductive
- Max. load: 30 kg
- Dimensions: 730 x 540 x 950/1350 mm
- Weight: 30 kg

## Hydraulic pump mayo instrument table

**Art.-no. 30001247**

Same as item 30001246, but with:

- Table top: 740 x 540 mm
- Dimensions: 780 x 540 x 950/1350 mm
- Weight: 31 kg



## Mayo instrument table

**Art.-no. 30001245**

- Mobile 4-leg construction made of stainless steel, seamlessly welded
- Stable, deep-drawn table top made of stainless steel, with rounded corners and edges, 640 x 440 mm
- Manual height adjustment from 830 - 1300 mm

- 4 double swivel castors Ø 75 mm, 2 of them electrically conductive
- Max. load: 30 kg
- Dimensions: 640 x 440 x 830/1300 mm
- Weight: 7 kg

# OR swivel chairs and OR stools



**Basic configuration: OR swivel chairs**

- Ergonomic, robust, hygienic and easy to clean
- Mobile five-castor construction made of stainless steel
- Infinitely variable height adjustment from 510 to 740 mm by foot control via chrome-plated safety foot bar
- Upholstery: Black, electrically conductive, with pull-out protection
- Round seat with ø 360 mm, thickness 100 mm / saddle seat with ø 420 mm, thickness 80 mm
- 5 double swivel castors ø 50 mm, 3 of them are electrically conductive
- With integrated foot brake in chassis
- Max. load: 140 kg

**OR swivel chair with round seat**

**Art.-no. 30001237**

Dimensions: 480 x 510 / 740 mm

Weight: 25 kg

**OR swivel chair with round seat and backrest**

**Art.-no. 30001238**

Backrest adjustable by 45 mm

Dimensions: 480 x 510 / 740 mm

Weight: 25 kg

**OR swivel chair with saddle seat**

**Art.-no. 30001239**

Dimensions: 480 x 510 / 740 mm

Weight: 30 kg

**OR swivel chair with saddle seat and backrest**

**Art.-no. 30001240**

Backrest adjustable by 45 mm

Dimensions: 480 x 510 / 740 mm

Weight: 30 kg



**Basic configuration: OR stools**

- Ergonomic, robust, hygienic and easy to clean
- Five-armed base made of stainless steel
- With precision spindle height adjustment from 510 to 710 mm
- Upholstery: Black, electrically conductive, with pull-out protection
- Round seat with ø 360 mm, thickness 70 mm
- Fixed version with 5 floor glides, 2 of them are electrically conductive
- Mobile version with 5 double swivel castors with interval braking ø 65 mm, 2 of them are electrically conductive
- Max. load: 140 kg

**OR stools with feet**

**Art.-no. 366050**

Dimensions: 600 x 510 / 710 mm

Weight: 7,6 kg

**Mobile OR stool**

**Art.-no. 366080**

Dimensions: 600 x 510 / 710 mm

Weight: 8 kg



# Instrument trolleys



**Basic configuration: mounted models**

- Mobile, jointless square tube frame construction made of stainless steel 20 x 20 mm
- Tube frame designed as a protective railing, top and bottom with removable stainless steel table tops
- Max. load: 20 kg per table top
- Drawer dimensions: 445 x 95 mm (only for variants with drawer)
- 4 abrasion-resistant, non contact-staining, corrosion-resistant antistatic double castors, ø 75 mm, 2 of them can be locked, with wall deflector
- Available in heights of 650, 800, 900 and 1,000 mm

**Basic configuration: unmounted models**

- Mobile, jointless square tube frame construction made of stainless steel 25 x 25 mm, cross braces 30 x 20 mm
- Tube frame designed as a protective railing, top and bottom with removable stainless steel table tops
- Max. load: 20 kg per table top
- 4 abrasion-resistant, non contact-staining, corrosion-resistant antistatic double castors, ø 75 mm, 2 of them can be locked



Description	Width x depth (in mm)	Height (in mm)	Art.-no.
Instrument trolleys (unmounted)	650 x 500	800	391075
	750 x 650	800	391275
	950 x 650	800	391475
Instrument trolleys (mounted)	675 x 475	650	30001353
		800	30001308
		900	30001309
		1000	30001310
	770 x 675	650	30001354
		800	30001311
		900	30001312
		1000	30001313
	990 x 675	800	30001314
		900	30001315
		1000	30001316
		1000	30001316
Instrument trolley with drawer (mounted)	675 x 475	650	30001355
		800	30001317
		900	30001318
		1000	30001319
	770 x 675	650	30001356
		800	30001320
		900	30001321
		1000	30001322
	990 x 675	800	30001323
		900	30001324
		1000	30001325
		1000	30001325
Instrument trolley with 3-sided stainless steel railing (mounted)	675 x 475	650	30001359
		800	30001326
		900	30001327
		1000	30001328

Description	Width x depth (in mm)	Height (in mm)	Art.-no.
Instrument trolley with 3-sided stainless steel railing (mounted)	770 x 675	650	30001360
		800	30001329
		900	30001330
		1000	30001331
	990 x 675	800	30001332
		900	30001333
		1000	30001334
Instrument trolley with drawer and 3-sided stainless steel railing (mounted)	675 x 475	650	30001357
		800	30001335
		900	30001336
		1000	30001337
	770 x 675	650	30001358
		800	30001338
		900	30001339
		1000	30001340
	990 x 675	800	30001341
		900	30001342
		1000	30001343
		1000	30001343
Instrument trolley with flush frame plate (mounted)	675 x 475	800	30001344
		900	30001345
		1000	30001346
	770 x 675	800	30001347
		900	30001348
		1000	30001349
	990 x 675	800	30001350
		900	30001351
		1000	30001352
		1000	30001352

# All-purpose and transport trolleys



**Basic configuration: All-purpose trolley**

- Mobile, jointless square tube frame construction made of stainless steel 20 x 20 mm
- Equipped with practical push handle
- Removable stainless steel table top with frame as a protective railing
- All trolleys are equipped with four robust wall deflectors
- Four swivel castors Ø 125 mm with thread guard, 2 lockable, electrically conductive
- Max. load: 20 kg per table top



**Basic configuration: Transport trolley**

- Stable, jointless round tube Ø 25 mm, made of stainless steel
- Fixed and welded shelves without joints and with an all round 10 mm raised profile edge, reinforced with 2 profiles on the underside
- All trolleys are equipped with rugged wall guards
- Four swivel castors Ø 125 mm with thread guard, 2 lockable, electrically conductive
- Max. load: 120 kg in total, 80 kg per shelf

Transport and multi-purpose trolleys must meet the highest standards when it comes to hygiene, economy, durability and ease of care. Therefore, all units are made of stainless steel due to the fact that it is even resistant to the most aggressive disinfectants, meaning that thorough cleaning is not an issue.

The jointless welding prevents residual dirt and germs from accumulating. All units have electrically-conductive castors.

The trolleys are suitable for universal use and particularly for transporting consumer goods up to 20 kg with a practical push handle.

**All-purpose trolley with two shelves**

**Art.-no. 30001361**

Dimensions (W x D x H):  
833 x 705 x 800 mm

**All-purpose trolley with two shelves**

**Art.-no. 30001362**

Dimensions (W x D x H):  
1053 x 705 x 800 mm

**All-purpose trolley with three shelves**

**Art.-no. 30001363**

Dimensions (W x D x H):  
1053 x 705 x 1020 mm

The transport trolley is the even more rugged answer to the all-purpose trolley. Loads of up to 120 kg (max. 80 kg per tray) can be transported without any problems.

The transport trolleys are suitable as serving trolleys and for transport of heavy equipment and objects.

**Transport trolley MSW 8 x 5-2**

**Art.-no. 30001367**

With two shelves 800 x 500 mm,  
dimensions (W x D x H):  
900 x 660 x 950 mm

**Transport trolley MSW 8 x 5-3**

**Art.-no. 30001368**

With three shelves 800 x 500 mm,  
dimensions (W x D x H):  
900 x 660 x 950 mm

**Transport trolley MSW 9 x 6-2**

**Art.-no. 30001369**

With two shelves 900 x 550 mm,  
dimensions (W x D x H):  
1000 x 650 x 950 mm

**Transport trolley MSW 9 x 6-3**

**Art.-no. 30001370**

With three shelves 900 x 550 mm,  
dimensions (W x D x H): 1000 x 650 x  
950 mm

**Transport trolley MSW 10 x 6-2**

**Art.-no. 30001371**

With two shelves 1000 x 600 mm,  
dimensions (W x D x H): 1100 x 700 x  
1010 mm

**Transport trolley MSW 10 x 6-3**

**Art.-no. 30001372**

With three shelves 1000 x 600 mm,  
dimensions (W x D x H): 1100 x 700 x  
1010 mm

# Hook systems



## Basic configuration: Hook trolley

- Mobile, jointless square tube frame construction made of stainless steel 25 x 25 mm
- Four antistatic swivel castors Ø 75 mm with wall deflectors, 2 of them lockable
- Vertical maximum load of hooks: 60 kg per hook
- Dimensions: 600 x 460 x 1438 mm
- Weight: 7 kg

Our comprehensive supply and disposal system consisting of hook trolley, wall-mounted hook rail and special accessories offers complete solutions in modular construction.

Assemble your hook system as you need it - with accessories for every conceivable task.

## Hook trolley, one side

**Art.-no. 30001261**

Models as described, with 3 pairs of hooks, for storage of up to three STE 1 baskets or comparable

## Hook trolley, both sides

**Art.-no. 30001262**

Models as described, with 3 pairs of hooks per side, for storage of up to six STE 1 baskets or comparable



## Basic configuration: Wall-mounted hook rail

- Square tube 30 x 20 mm with welded hooks and two synthetic caps
- Hooks with vertical distance of 350 mm
- Vertical maximum load of hooks: 60 kg per hook
- Package unit: piece (for attachment of accessories one pair is needed)

## Wall-mounted hook rail with 2 hooks

**Art.-no. 30001263**

Models as described, with two hooks, length: 700 mm, weight: 1 kg

## Wall-mounted hook rail with 3 hooks

**Art.-no. 30001264**

Models as described, with three hooks, length: 1050 mm, weight: 1,4 kg

## Wall-mounted hook rail with 4 hooks

**Art.-no. 30001265**

Models as described, with four hooks, length: 1400 mm, weight: 1,4 kg

## Wall-mounted hook rail with 5 hooks

**Art.-no. 30001266**

Models as described, with five hooks, length: 1750 mm, weight: 2,2 kg

## Wall-mounted hook rail with 7 hooks

**Art.-no. 30001267**

Models as described, with seven hooks (vertical distance 240 mm) for attachment of seven STE 1/2 baskets, length: 1600 mm, weight: 2,5 kg



# Accessories for hook systems

## Sterile equipment baskets

Can be suspended in wall-mounted hook rail and hook trolley. For transport, storage and sterilisation. Made of electropolished stainless steel. Max. Load: 20 kg Weight: 2 kg

### Sterile equipment basket 1/2 STE

**Art.-no. 30001268**  
Model as described, dimensions (W x D x H): 575 x 280 x 135 mm

### Sterile equipment basket 1 STE

**Art.-no. 30001269**  
Model as described, dimensions (W x D x H): 575 x 280 x 265 mm

### Sterile equipment basket 1 STE with handle opening

**Art.-no. 30001270**  
Model as described, but with handle opening at long side 300 x 120 mm, dimensions (W x D x H): 575 x 280 x 265 mm

### Sterile equipment basket 2 STE

**Art.-no. 30001271**  
Model as described, but with maximum load of 30 kg, dimensions (W x D x H): 575 x 560 x 265 mm



## Universal holder for synthetic and wire containers

For holding of containers with size 3, 4 and 5, made of electropolished stainless steel. Dimensions (W x D): 550 x 37 mm Weight: 1 kg

### For synthetic and wire container 3

**Art.-no. 30001272**  
Model as described, for holding of 3 containers with size 3, height: 247 mm

### For synthetic and wire container 4

**Art.-no. 30001273**  
Model as described, for holding of 4 containers with size 4, height: 170 mm

### For synthetic and wire container 5

**Art.-no. 30001274**  
Model as described, for holding of 5 containers with size 5, height: 124 mm



## Synthetic container

Can be suspended in universal holders for containers, made of polypropylene, stackable.

### Synthetic container size 3, blue


**Art.-no. 30001294**  
Dimensions: 210 x 350 x 200 mm

### Synthetic container size 4, blue

**Art.-no. 30001295**  
Dimensions: 150 x 230 x 125 mm

### Synthetic container size 5, blue

**Art.-no. 30001296**  
Dimensions: 100 x 170 x 75 mm



## Wire container

Can be suspended in universal holders for containers, made of stainless steel, stackable.

### Wire container size 3

**Art.-no. 30001297**  
Dimensions: 210 x 350 x 200 mm

### Wire container size 4

**Art.-no. 30001298**  
Dimensions: 150 x 230 x 125 mm

### Wire container size 5

**Art.-no. 30001299**  
Dimensions: 100 x 170 x 75 mm



# Accessories for hook systems



**Multi-purpose tray**

**Art.-no. 30001275**


- Holds folded paper
- Jointless stainless steel finish
- Dimensions: 610 x 326 x 260 mm
- Weight: 1,5 kg



**Shoe rack**

**Art.-no. 30001276**

- Holds a maximum of 6 pairs of shoes
- Jointless finish made of electropolished stainless steel
- With 12 hanging bars
- Closed wire loops with wire thickness 6 mm
- Dimensions: 600 x 263 x 180 mm
- Weight: 1 kg



**Hose and abdominal cloth rack**

**Art.-no. 30001277**

- Holds anaesthesia hoses and abdominal cloths
- Jointless stainless steel finish
- With 11 welded hooks, arranged offset
- Dimensions: 585 x 94 x 180 mm
- Weight: 1,5 kg



**Multi-purpose equipment holder**

**Art.-no. 30001278**

- For example for anaesthesia units, operating theatre accessories
- Jointless stainless steel finish, slide rail produced of solid material 25 x 10 mm
- Dimensions: 600 x 48 x 180 mm
- Weight: 1,5 kg



**Stainless steel housing for acrylic glass chutes (with 3 chutes)**

**Art.-no. 30001279**

- Stackable chutes
- Dimensions: 601 x 172 x 245 mm
- Weight: 5 kg

**Stainless steel housing for acrylic glass chutes (with 5 chutes)**

**Art.-no. 30001280**

- Stackable chutes
- Dimensions: 601 x 172 x 171 mm
- Weight: 4,5 kg



**Holder for stainless steel container**

**Art.-no. 30001281**

- Jointless stainless steel finish
- With storage bar for lids

**Stainless steel container (100 mm)**

**Art.-no. 30001282**

- Dimensions: 530 x 325 x 100 mm
- Weight: 1,4 kg

**Stainless steel container (200 mm)**

**Art.-no. 30001283**


Same as item 30001282, but with dimensions: 530 x 325 x 200 mm, weight: 1,4 kg

**Lid for stainless steel container**

**Art.-no. 30001284**

- Dimensions: 530 x 325 mm
- Weight: 1,1 kg


# Accessories for hook systems



**Syringe and cannula tray**

**Art.-no. 30001285**

- Equipped with 4 acryl glass subdivisions, each can be adjusted variably
- Jointless stainless steel finish
- Dimensions: 580 x 350 x 275 mm
- Weight: 9 kg



**Suture material tray**

**Art.-no. 30001286**

- With 4 shelves, holds a maximum of 20 commercially available small packages
- Jointless stainless steel finish
- Dimensions: 600 x 173 x 290 mm
- Weight: 6 kg



**Glove dispenser**

**Art.-no. 30001287**


- Holds a maximum of 50 pairs of 3 glove sizes
- Jointless stainless steel finish
- Dimensions: 580 x 320 x 260 mm
- Weight: 9 kg



**Writing surface**

**Art.-no. 30001288**

- As writing support for nursing staff
- Jointless stainless steel finish
- Writing surface with front bevel and middle shelf
- Dimensions: 580 x 396 x 244 mm
- Weight: 6 kg



**Storage for urinal bottles**

**Art.-no. 30001289**

- Holds 5 urinal bottles
- Jointless stainless steel finish
- Dimensions: 580 x 276 x 230 mm
- Weight: 4 kg

**Storage for 5 bedpans with lids**

**Art.-no. 30001290**

- Jointless stainless steel finish
- Dimensions: 675 x 420 x 230 mm
- Weight: 6 kg



**Support surface with bracket (350 x 205)**

**Art.-no. 30001291**

- With permanently-welded tray made of stainless steel
- Jointless stainless steel finish
- Load capacity with evenly distributed load: 50 kg
- Dimensions: 695 x 350 x 205 mm
- Weight: 5 kg

**Support surface with bracket (450 x 240)**

**Art.-no. 30001292**

Same as item 30001291, but with

- Dimensions: 695 x 450 x 240 mm
- Weight: 6 kg

**Support surface with bracket (600 x 275)**

**Art.-no. 30001293**

Same as item 30001291, but with

- Dimensions: 695 x 600 x 275 mm
- Weight: 8 kg

# Accessories for hook systems



**Wire mesh basket with folding handle (240 mm)**

**Art.-no. 30001300**

- Wire mesh basket deep-drawn, closed base frame with crossbars for high loading capacity
- 2 folding handles inside
- Jointless electropolished stainless steel finish
- Handle bar: 130 x 40 mm
- Dimensions: 240 x 250 x 50 mm

**Wire mesh basket with folding handle (480 mm)**

**Art.-no. 30001301**

Same as item 30001300, but with dimensions: 480 x 250 x 50 mm



**Strainer**

**Art.-no. 30001302**

- Jointless electropolished stainless steel finish
- Closed base frame, stackable
- Recess lengthwise in upper frame serves as carrying bar
- 4 U-shaped feet welded to the base frame
- Fits in polypropylene containers
- Dimensions: 240 x 250 x 50 mm

**Polypropylene container with closed lid**

**Art.-no. 30001303**

- To take up strainer
- Can be sterilised up to +134 °C
- Dimensions: 300 x 225 x 125 mm

**Polypropylene container with slotted lid**

**Art.-no. 30001304**

Same as item 30001303, but with slotted lid



**Pallet trolley**

**Art.-no. 30001307**

- Mobile, jointless square tube frame construction, 25 x 25 mm, made of stainless steel
- Holds two 1 STE baskets or one 2 STE basket
- With moulded push handle
- 4 steering castors Ø 125 mm with double ball bearing, electrically conductive, with thread guard and wall guard, 2 of them with brakes
- Max. load: 100 kg
- Dimensions: 675 x 815 x 900 mm

# Further OR furniture



**Bowl stand with heating**

**Art.-no. 30001243**

- Bowl: made of stainless steel, with positioning ring and push handle, volume 5,5 litres
- Heating: radiator 230 V, 50/60 Hz 600 W, infinitely variable from 30-60 °C with lock catch at 40 °C, temperature setting in 10 °C increments, with LED indicator, control accuracy +/- 2 °C, 5 m connection cable and Schuko plug
- Frame: mobile, 5-foot design with tip safety made of stainless steel, height can be adjusted with one hand from 1004 to 1275 mm
- 5 double castors Ø 75 mm, 3 of them electrically conductive, 2 of them lockable
- Dimensions: 630 x 630 x 1004 / 1275 mm



**Bowl stand, single**

**Art.-no. 30001241**

- Made of stainless steel
- With one bowl (6 litres)
- 5 double castors Ø 75 mm, 3 of them electrically conductive, 2 of them lockable
- Dimensions: 635 x 670 x 865 mm

**Bowl stand, double**

**Art.-no. 30001242**

Same as item 30001241, but with two bowls and dimensions 635 x 700 x 865 mm



**Infusion stand**

**Art.-no. 30001232**

- Made of stainless steel
- Height-adjustable from 1500 to 2500 mm using screw
- Incl. stainless steel bottle holder with 4 safety infusion hooks and drip vessel and holder
- Max. load per hook: 2 kg
- 5 double castors Ø 75 mm, 2 of them lockable
- Dimensions: 630 x 630 x 1500/2500 mm

**Infusion stand, one-hand model**

**Art.-no. 30001252**

Same as 30001232 but with one-hand height adjustment from 1500 to 2350 mm and dimensions 630 x 630 x 1500/2350 mm



## Bucket stand

**Art.-no. 364010**

- With bucket 15 litres
- Made of stainless steel
- 4 double castors  $\varnothing$  65 mm, two of them electrically conductive
- With bumper ring
- Dimensions ( $\varnothing$  x H): 420 x 400 mm

## Bucket stand with push handle

**Art.-no. 364020**

Same as item 364010, but with

- Push handle
- Dimensions ( $\varnothing$  x H): 420 x 400 mm (push handle H= 860 mm)

## Lid for bucket stand

**Art.-no. 364030**

- For stainless steel bucket (15 litres)



## OR step stool, single, height 130 mm

**Art.-no. 30001233**

- Stable 4-foot design without joints made of stainless steel
- With removable, non-slip, corrugated step surface
- 4 electrically conductive, non-discolouring rubber feet
- Max. Load: 125 kg
- Dimensions: 600 x 350 x 130 mm

## OR step stool, single, height 220 mm

**Art.-no. 30001234**

Same as item 30001233, but with h= 220 mm

## OR step stool, single, height 300 mm

**Art.-no. 30001235**

Same as item 30001233, but with h= 300 mm

## OR step stool, double

**Art.-no. 30001236**

Same as item 30001233, but with two steps, dimensions 600 x 350 x 220/440 mm



## Transport trolley, single

**Art.-no. 30001249**

- Mobile, jointless tubular frame design made of stainless steel, with hanging fixture which holds commercially-available laundry or waste sacks with 60/70 l volume
- With moulded push handle
- With 4 double castors  $\varnothing$  75 mm, two of them lockable and electrically conductive
- Dimensions: 575 x 350 x 777 mm

## Transport trolley, double

**Art.-no. 30001250**

Same as item 30001249, but for two laundry or waste sacks, dimensions 1020 x 350 x 777 mm

## Synthetic lid for waste trolley

**Art.-no. 30001251**

- Dimensions: 403 x 360 mm



## Hand and arm surgical table

**Art.-no. 30001244**

- Mobile 4-foot design made of stainless steel and welded without joints
- Manually height adjustable from 650 to 950 mm
- Table top with radiolucent upholstery, height 30 mm, can be removed or tilted to all sides with ball joint and can be tightly fixed in place with clamping lever
- Max. load: 15 kg
- 4 abrasion-resistant, non contact-staining, corrosion-resistant double castors,  $\varnothing$  75 mm, two of them electrically conductive
- Dimensions: 450 x 450 x 650 / 950 mm



## OR shoe rack, mobile

**Art.-no. 30001253**

- Mobile 4-foot design made of stainless steel and welded without joints
- With 6 stainless steel tube frame levels as support for 24 pairs of operating theatre shoes of all sizes
- 4 abrasion-resistant, non contact-staining, corrosion-resistant double castors,  $\varnothing$  75 mm, two of them lockable
- Dimensions: 970 x 580 x 1600 mm

## OR shoe rack, wall-mounted

**Art.-no. 30001254**

- Stainless steel frame construction, welded without joints
- With 7 stainless steel tube frame levels as support for 30 pairs of operating theatre shoes of all sizes
- Dimensions: 894 x 345 x 1450 mm





**Basic configuration: Examination couch**

- Stable 4-foot design made of powder coated steel sheet 30 x 30 mm
- Adjustable head section up to +35°
- Upholstery: Artificial leather covered foam, height = 50 mm
- With spacers between upholstery and frame
- Incl. paper roll holder stick
- With four synthetic floor glides
- Dimensions: 1950 x 650 x 650 mm
- Max. load: 225 kg



**Basic configuration: Examination couch made of stainless steel**

- Stable 4-foot design made of stainless steel 30 x 30 mm
- Adjustable head section up to +35°
- Upholstery: Artificial leather covered foam, height = 70 mm
- With spacers between upholstery and frame
- Incl. paper roll holder stick
- With four synthetic floor glides
- Dimensions: 1950 x 650 x 650 mm
- Max. load: 225 kg

**Examination couch**

**Art.-no. 111419**

Model as described

**Examination couch**

**Art.-no. 111519**

Model as described but with upholstery height: 70 mm

**Mobile examination couch**

**Art.-no. 111519F**

Model as described, but with:

- Upholstery height: 70 mm
- Longitudinal and transversal bracing
- 4 steering castors Ø 100 mm, 2 of them lockable, 1 of them direction-adjustable, 1 of them free-wheeling
- 4 wall guards

**Examination couch made of stainless steel**

**Art.-no. 171519**

Model as described

**Mobile examination couch made of stainless steel**

**Art.-no. 171519F**

Model as described, but with:

- Longitudinal and transversal bracing
- 4 steering castors Ø 100 mm, 2 of them lockable, 1 of them direction-adjustable, 1 of them free-wheeling
- 4 wall guards



**Basic configuration: Mobile stools standard**

- Polished five-armed aluminium base
- Gas spring supported height adjustment from 420 to 550 mm
- Five double castors with interval braking  $\varnothing$  50 mm
- Upholstery: Artificial leather covered foam according to colour chart
- Seat:  $\varnothing$  360 mm, h = 60 mm
- Base:  $\varnothing$  600 mm



**Basic configuration: Mobile stools with arm/back rest**

- Incl. height adjustable arm or back rest
- Polished five-armed aluminium base
- Gas spring supported height adjustment from 420 to 550 mm via release lever under seat
- Five double castors with interval braking  $\varnothing$  50 mm
- Upholstery: Artificial leather covered foam according to colour chart
- Seat:  $\varnothing$  360 mm, h = 60 mm
- Base:  $\varnothing$  600 mm

**Mobile stool**

**Art.-no. 206020**

Model as described, height adjustment via release lever under seat

**Mobile stool**

**Art.-no. 206030**

Same as 206020, but with height adjustment from 510 to 690 mm

**Mobile stool with saddle seat**

**Art.-no. 206040**

Same as 206020, but with saddle seat  $\varnothing$  400 mm

**Mobile stool with saddle seat**

**Art.-no. 206050**

Same as 206040, but with height adjustment from 510 to 690 mm

**Mobile stool with foot release**

**Art.-no. 206060**

Model as described, height adjustment from 435 to 565 mm via foot release under base, with double castors with  $\varnothing$  65 mm

**Mobile stool with foot release**

**Art.-no. 206070**

Same as 206060, but with height adjustment from 525 to 705 mm

**Mobile stool with back rest**

**Art.-no. 206024**

Model as described

**Mobile stool with back rest**

**Art.-no. 206034**

Same as 206024, but with height adjustment from 510 to 690 mm

**Mobile stool with back rest and base ring**

**Art.-no. 206035**

Model as described, but with base ring and height adjustment from 510 to 690 mm

**Electric conductivity**

**Art.-no. 206012**

Surcharge for electrically conductive castors and covering for stool models 206020, 206030, 206060 & 206070

**Electric conductivity**

**Art.-no. 206013**

Surcharge for electrically conductive castors and covering for stool models 206024, 206034 & 206035

**Electric conductivity**

**Art.-no. 206014**



















Surcharge for electrically conductive castors and covering for stool models 206040 & 206050

# Upholstery colours

**For examination couches and mobile stools**

- Ergonomic implementation of upholstery parts
- Rounded rough edges
- High quality artificial leather
- Seamlessly manufactured
- Optimal possibilities to clean and to disinfect
- Various colourings without any extra charge



 -50	 -56	 -63	 -67	 -75
White	Beige	Lime	Lemonade	Agave
 -73	 -61	 -65	 -53	 -58
Apricot	Carnelian	Cherry	Pearl	Stone
 -55	 -62	 -74	 -80	 -81
Dolphin	Light Blue	Atoll	Royal	Purple
 -107	 -83	 -96		
Graphite grey	Black	Black, electrically conductive		

# Further products for the OR

**RooSy.** The modular Room System of medifa.

The prefabricated modular room system from medifa is specially designed for the high demands in the operating room.

The perfectly coordinated wall, door and ceiling system guarantees the highest standards with maximum flexibility. medifa undertakes the complete planning, production and installation in-house and sees itself as a complete supplier of turnkey operating theatres.

Learn more:



# Further products for the OR

## Mobile operating tables of medifa.

Functionality and comfort for more efficiency in the operating room.

The mobile operating tables from medifa enable optimal, simple and safe patient positioning in the operating theatre.

The operating tables are suitable for all surgical disciplines, requirements and applications in clinics, hospitals and outpatient operating centres.

Learn more:



# medifa



The medifa healthcare group unites competence concerning operating rooms under one roof.

medifa is an owner-run and internationally active complete solution provider in the medical technology.

The product and service portfolio covers all operating room requirements from modular wall, door, ceiling and cabinet systems (RooSy® Room Systems) for the turnkey operating room, as well as mobile operating tables, operating table accessories and operating lights, right up to examination chairs and medical

furniture. In addition, one of the largest German OEM suppliers in the medical technology sector also belongs to medifa.

medifa GmbH & Co. KG (Finnentrop), medifa hygienic rooms GmbH (Ötigheim) and medifa metall und medizintechnik GmbH in Rastatt are all part of this group of companies. The medifa healthcare group employs approximately 350 employees at three locations within Germany. The planning, engineering and production


take place exclusively in Germany, and with over 90% added value in the group of companies, the group guarantees fast reaction times.

The slogan “we care” means for us that we consistently orientate ourselves according to the requirements of our customers. Customers all over the world can rely on top-class consultation, processing and quality from a single source, made in Germany according to certified national and international quality standards.

## Products

- Modular operating room systems: walls, doors, ceilings and cabinets (RooSy)
- Mobile operating tables and accessories
- Operating lights
- Examination and treatment chairs
- Medical furniture and trolleys
- Customized medical solutions

## medifa worldwide



we care.





Version 01-20-EN  
Data subject to change without notice.

[www.medifa.com](http://www.medifa.com)



**medifa GmbH & Co. KG**

Industriestraße 5  
57413 Finnentrop  
Germany  
fon +49 2721 71 77 -0  
fax +49 2721 71 77 -255  
[info@medifa.com](mailto:info@medifa.com)

**we care.**